MATERIALS AND STRUCTURES SYMPOSIUM (C2) Interactive Presentations (IP)

Author: Mr. Ruven Spannagel DLR, German Aerospace Center, Germany, ruven.spannagel@dlr.de

Ms. Ines Hamann
DLR (German Aerospace Center), Germany, Ines.Hamann@dlr.de
Dr. Jose Sanjuan
DLR, German Aerospace Center, Germany, jose.sanjuan@dlr.de
Dr. Felipe Guzman
DLR (German Aerospace Center), Germany, Felipe.Guzman@dlr.de
Prof. Claus Braxmaier
ZARM University of Bremen, Germany, Claus.Braxmaier@zarm.uni-bremen.de

DIMENSIONAL STABILITY INVESTIGATION OF LOW CTE MATERIALS AT TEMPERATURES FROM 140 K TO 250 K USING A HETERODYNE INTERFEROMETER

Abstract

Light weight materials with excellent dimensional stability are increasingly needed in space based applications such as telescopes, optical benches, and optical resonators. Glass-ceramics and composite materials can be tuned to reach very low coefficient of thermal expansion (CTE) at certain temperatures, including room temperature and cryogenics, where a growing number of instruments in scientific and earth observation space missions are operated. Very accurate setups are needed to determine the CTE of such materials. With our laser-interferometric dilatometer setup we are able to measure CTEs of a large variety of materials in the temperature range of 140 K to 250 K. Special mirror mounts with a thermally compensating design enable measurements of the expansion of cylindrical tube-shaped samples using a heterodyne interferometer with demonstrated noise levels in the order of 10 pm/ $\sqrt{\rm Hz}$. The temperature variation of the sample is obtained by a two stage controlled heating/cooling setup where a pulse tube cooler and electric heaters apply small amplitude temperature signals to cool/heat the sample radiatively in order to reach a homogeneous temperature over the whole sample. A carbon fiber reinforced polymer (CFRP) sample was selected to run CTE measurements, achieving results in the 10^{-8} K⁻¹ range including all known uncertainties. The limitations of our setup have been identified and the largest uncertainty contribution has been determined to be tilt-to-length coupling of the sample due to temperature variations. Several improvements are currently underway to minimize our uncertainty budget. New results with the enhanced setup will be presented